

## *Basket Sinker Specifications*

### Hand Fabricated Style



Catalog Number	Max Capsule Size	Inside Dia.	Min. Inside Length	Weight
BSK008-01	000	.62 inch	1.06 inches	5.0 gm
BSK008-JP	1	.51 inch	.90 inches	3.6 gm

Lid ----- 316 Stainless Steel  
 Mesh ----- 316 Stainless Steel, 8 Mesh, Wire Cloth, .028 Wire Dia.  
 Finish ----- Electropolish

#### **Ordering Information:**

QLA Basket Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your Quality Lab Accessories Sales Representative to order others sizes we carry, or refer to our web site at [www.qla-llc.com](http://www.qla-llc.com) to see all our product line.