




## *Capsule Sinker Specifications Helix Coil Type*



### CAPWHT-4S

Material	316 L Stainless Steel
Finish	Electropolished
Inside Length	.840 inches
Inside Diameter	.385 inches
Wire Diameter	.047 inches
Pitch	.280 inches
End Coil Inside Diameter	.190 inches
Maximum Width	.60 inches
Number of Coils	5
Maximum Capsule Size	0

### **Ordering Information:**

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your Quality Lab Accessories Sales Representative to order others sizes we carry, Or refer to our web site at [www.qla-llc.com](http://www.qla-llc.com) to see all our product line.