



Capsule Sinker Specifications

USP Compliant Alternative Capsule Sinker with *Magnetic Qualities*



CUSBSK-JPMAG (For Automated Retrieval)

Material (Magnetic)	316L Stainless Steel
Overall Length	25mm – 26mm
Coil Spacing (Pitch)	3.0mm – 3.5mm
Clasp Lid & Bottom Spacing	3.5mm – 4.00mm
Inside Diameter	12.0mm +/- 0.2mm
Acid Resistant Wire Clasp	316L Stainless Steel
Acid Resistant Wire Support (10ea)	316L Stainless Steel

Ordering Information:

QLA Capsule Sinkers come in a wide range of sizes to fit all dissolution applications, contact your Quality Lab Accessories sales representative to order others sizes we carry, or refer to our web site at www.qla-llc.com to see all our product line.